Redyrob® Trans Plus

en	Closed, controllable wound drainage system	2
bg	Затворена, регулируема система за дрениране на рани	4
cs	Uzavřený regulovatelný systém pro drenáž rány	6
da	Lukket, regulerbart sådrænesystem	8
de	Geschlossenes, kontrollierbares Wunddrainagesystem	10
el	Κλειστό, ελεγχόμενο σύστημα παροχέτευσης τραύματος	12
es	Sistema de drenaje de heridas cerrado y regulable	14
et	Suletud, kontrollitav haavadrenaaži süsteem	16
fi	Suljettu, kontrolloitava haavandreneerausjärjestelmä	18
fr	Système de drainage des plaies fermé et contrôlable	20
hr	Zatvoren, kontrolirani sustav za drenažu rane	22
hu	Zárt, kontrollálható sebdrenázs rendszer	24
id	Sistem drainase luka tertutup dan terkendali	26
it	Sistema di drenaggio della ferita controllabile e chiuso	28
lt	Uždara, kontroliuojama žaizdų drenavimo sistema	30
lv	Slēgta, regulējama brūces drenāžas sistēma	32
nl	Gesloten regelbaar wonddrainagesysteem	34
no	Lukket, kontrollerbart sådreneringssystem	36
pl	Zamknięty, regulowany system drenażu rany	38
pt	Sistema de drenagem de feridas controlável fechado	40
ro	Sistem închis de drenaj pentru plăgi, cu aspirație controlată	42
sk	Uzatvorený kontrolovateľný drenážny systém na rany	44
sl	Zaprto, nadzorovan sistem za drenažo rane	46
sv	Slutet, kontrollerbart sådräneringssystem	48
tr	Kapalı, kontrol edilebilir yara drenaj sistemi	50

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


en - Redyrob® TransPlus - 600 ml

Device description:

Closed adjustable wound drainage system for connection to wound drainage

– 600 ml secretion container, not made with PVC

– Vacuum regulator

– Quantitative vacuum gauge

The closed circuit and electively pressure adjustable wound drainage system Redyrob® eliminates the situation to have to substitute containers during the drainage process and therefore protects patients from ascending infective organisms and the staff from inadvertently disseminated infectious material.

Intended use:

Wound drainage systems are used for the postoperative drainage of blood and secretions from wounds and body cavities for hematoma and seroma prophylaxis.

Indication:

Wound drainage systems are used for the postoperative drainage of blood and secretions from wounds and body cavities for hematoma and seroma prophylaxis. They are used for initial wound healing as part of holistic postoperative wound care.

Redyrob® wound drainage systems are particularly indicated for orthopedic procedures such as hip and knee joint operations.

By draining the wound secretion, the wound cavity is reduced and the wound surfaces are drawn together. Vacuum drainage thus ensures better adaptation and stabilization of the wound surfaces and reduces the risk of infection. Wound healing is accelerated.

Contraindication:

In general, the contraindications to wound drainage under vacuum must be taken into account.

- Severe bleeding may occur if the drainage system is inserted close to open cancellous bone
- Vacuum drains must not be placed in direct contact with the brain, liver, spleen or other parenchymatous organs
- There is a risk of severe blood loss when inserting vacuum drains into trabecular bone wounds
- Not to be used in patients with physiological blood coagulation disorders

Intended users:

A medical degree must always be named as the person responsible, as the placement of the drainage is reserved for this person. If necessary, further manipulation/application is then carried out by qualified nursing staff.

Patient population:

Wound drainage systems are part of the basic surgical and clinical equipment and are used for (emergency) patients without restriction based on gender or age. A potential physiological limitation is a blood clotting disorder.

Apart from this physiological contraindication, wound drainage systems are not subject to any specific assignment to a patient population.

Adverse side effects and complications:

The following risks are known from years of use of the products and are considered as part of risk management:

- In the wound area: redness, swelling, production of purulent secretion
- Ingrowth of the suction drain after a longer period of stay and subsequent unpleasant removal for the patient
- General therapy-related side effects: hematoma/seroma formation, bleeding, discomfort, pain, infection
- Iatrogenic trauma due to misplacement or displacement of the drain
- Obstruction of the drain
- Increased rate of infection due to movement of bacteria across the drain

Warnings:

- Use exclusively for the intended purpose
- Do not use if packaging is damaged
- Disposable product, do not empty the contents and do not reuse the systems under any circumstances. Reuse of single-use products poses a potential risk to the patient or user. It may lead to contamination and/or impairment of functionality. Contamination and/or impaired functionality of the product may result in patient injury, illness or death
- Use an adequate number of wound drains to ensure that the wound or body cavity is draining properly
- Ensure correct placement of the drain(s) in the wound area
- Adequate securing of the drain to avoid disconnection or displacement of the system
- Ensuring that the wound is hermetically sealed
- Check the connection of the drain and drain connector before use
- Checking the operational readiness of pre-evacuated wound drainage systems: Checking the vacuum display using the quantitative vacuum gauge
- The type and amount of wound secretion collected should be checked regularly
- In case of acute bleeding: stop treatment immediately
- Regular checking of the proper functioning of the drainage
- Replace the wound drainage container as soon as there is no more vacuum (vacuum gauge display shows "low")
- Danger! Risk of stabbing injuries!
- Follow national needlestick guidelines and regulations
- Observe national reporting requirements in the event of incidents. Any serious incident should be reported to the manufacturer and the local regulatory authority.
- Dispose of used products in a contamination-free manner in accordance with national disposal guidelines! Infection risk!

Directions for use:

1. Check that quantitative vacuum gauge shows vacuum in vessel.
2. Remove the red Combi-Stopper from Luer Lock of the Redyrob® TransPlus.

CAUTION: Do not connect anything to the red Combi-Stopper! The Combi-Stopper prevents wound secretion flow!

Connect Redyrob® TransPlus to the wound drainage using the LuerLock connector.

Set the desired vacuum according to the following table on the Redyrob® TransPlus.

Functions of control positions 0–3:

Position 0: The inlet is closed. Thus, wound secretion cannot flow away.

This may cause tamponade of the wound or the formation of hematomas

Position 1: Low vacuum is applied to the wound

Position 2: Medium vacuum is applied to the wound

Position 3: High vacuum is applied to the wound

In positions 1–3 wound secretion is sucked into the Redyrob® TransPlus vessel.

3. Turn the cap to the desired vacuum level. The system is active.

If the vacuum gauge shows "low" and vacuum drainage is still needed, connect another Redyrob® TransPlus.

Duration of use

General statement not possible as drainage duration is guided by clinical observation.

If the suction bottle is full or contains air, vacuum drainage is no longer possible without the addition of a Redyrob® TransPlus system.

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Boston Scientific (Australia) Pty Ltd
MEDINORM Medizintechnik GmbH