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
ENGLISH

LVS: Low Vacuum Wound Drainage System - 600 mL

Application

Sterile, disposable wound drainage system to be used for the removal of sero-sanguineous fluid from the wound site in order to assure proper wound approximation and healing, to alleviate dead space, and to reduce potential complications.

Precautions

- introduce a sufficient number of wound drains to ensure adequate drainage
- ensure the wound drain(s) is (are) placed correctly
- ensure an airtight wound seal
- ensure the vacuum indicator (green bellow (4)) is in the "max" position
- the nature and amount of collected wound secretion should be checked (2) on a regular basis and indicated on the bottle (3)
- pronounced haemorrhaging may occur if the drainage system is inserted next to open cancellous bone
- avoid direct contact between the wound drain and the brain, or the intestinal, thoracic, as well as any other parenchymal organ
- STERILE; unless package opened or damaged

Storage

- to be stored in its packaging at a temperature between 0°C (32°F) and 25°C (77°F), with a relative humidity between 50% and 60%
- peak temperatures of 35°C (95°F) should never be exceeded
- prevent placing the product in direct sunlight

Directions for use

1. use aseptic techniques when opening sterile package
2. after removing the needle cover (19) (caution: danger of injury), introduce the wound drain (17) with the needle (18) in the wound area going from the inside out, until the black marker line (16) is at skin level
3. cut the drain -close to the needle- at a 45° angle for easy introduction in the single- or Y-drain connector (14)
4. for single-drain connectors, go to step 7
5. a middle perforated drain can be cut in 2 parts
6. unscrew the Medifix drain connector (15) from the needle and connect the second part of the drain to the needle, before repeating steps 2 and 3
7. cut the single- or Y-drain connector at the appropriate drain size indicator (14-a or 14-b)
8. introduce the wound drain into the single- or Y-drain connector, check the correct fixation
9. squeeze the pressure regulator (7) using index finger and thumb
10. connect the connecting tube (11) with the Luer-Lock (8) on the bottle (1) (caution: loose connections can lead to loss of vacuum)
11. remove the security-clip (5) above the "off/on" lever (6)
12. open the lever from "off" to "on", pushing it slightly upwards to the left and then to the right; after releasing the pressure regulator, it will start functioning

THE SYSTEM IS NOW IN USE and can be attached to the bed, in any position, by using the bed strap (10).

Connection of a replacement bottle

1. if a replacement is necessary, e.g. when the vacuum indicator is in the "minimum" position, close the sliding clamp (12) on the connecting tube
2. close the lever by moving it from the "on" to the "off" position
3. disconnect the connecting tube by unscrewing the Luer-Lock; backflow into the patient will be prevented by the anti-reflux valve (13) (if applicable)
4. close the full bottle by using the screw-cap (9)
5. discard the full system, properly sealed, according to normal hospital policy
6. take a replacement bottle and repeat the steps from "directions for use", point 9 - 12
7. open the sliding clamp on the connecting tube

Disconnection

1. after the surgeon (or any other authorized person) has taken the decision to stop the wound drainage, close the sliding clamp on the connecting tube and close the lever by moving it from "on" to the "off" position
2. before removing the drain from the patient's side, ensure the pressure difference inside and outside the patient will be eliminated, by letting the wound relax for 10 - 20 minutes
3. repeat steps 3, 4, and 5 of "connection of a replacement bottle"

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